

A full-featured, easy-to-use communication app for the iPad.

Tips

From author SLP Gisele Morin-Connolly



Success with Core Word Pages

- Designed for students with complex communication needs
- As the core vocabulary expands, the general layout remains consistent
- Ideally students with fine motor ability start with a 36-button layout with only a few messages visible as you teach language
- As student progresses make more messages visible begin at level student can physically access messages (4, 9, 16, or 25)
- Don't expect them to begin "talking" immediately
- Allow the student "wait time" as they may need more time to process and respond
- Take every opportunity to model with the same layout using your own iPad or the student's iPad



Turn your iPad into a GoTalk!

GoTalk[®] NOW with Core Words GoTalk[®] NOW Plus with Core Words

To order call 1-800-327-4269

.\$99.99 .\$169.99

FEATURES

- Easy-to-use page editor allows you to delete, crop, and scale images quickly.
- Use images from your iPad camera, photo library, or built-in Internet image search.
- 1-36 message locations per page plus up to 8 core messages!
- Randomize a message location i.e. Hi, Hello, Hey, or Hola.
- Add music or video for engagement, step-by-step video modeling, or social skills.
- Record your own speech or use free built-in iOS voices.
- Choose an Acapela high-quality textto-speech voice, starting at just 99¢ each (20+ languages).

www.AttainmentCompany.com

Available on the App Store

Switch Users

Auditory cues make GoTalk[®] NOW perfect for switch

alk	greetings		calendar	
? questions	describing	Lunch	about me	
Y ∦ need	math	iliteracy	social studies	
categories	art	science	Mar Ar music	
🟠 🕽 🔒 Home 🏂 👀				

users. Customize your scan settings for use with 1 or 2 switches. It's the only communication app to provide private auditory cues.

3 styles of personal communication pages



Standard

With standard pages, press any communication location for speech, music, video, or to jump to another page.



Express

Express pages build recorded or text-tospeech messages in a speech bar at the top or bottom of the screen. Press the bar to play the messages together in a sequence.



Scene

Build scene pages around a single photo or a pair of images. Arrange visible or invisible hotspots over people or objects in the photo, and the hotspots play speech, music, or videos.



Offers the complete **GoTalk**[®] **NOW** App without downloading any in-app purchases. Incorporates all the features of **GoTalk**[®] **NOW**.



())

SymbolStix library

Ready-Set-Communicate

Acapela[®] text-to-speech



www.AttainmentCompany.com 1-800-327-4269